

**A COMMON FACTOR FOR THE SAMPLE THAT AFFECTS THE EXAMINATION'S PERFORMANCE IN DISEASE SECTION**

UNIT	TEST	SAMPLE	FACTOR	CONSEQUENCES
Tibi/Kusta	MTB Identification and Drug susceptibility testing (DST)	Sputum	<ul style="list-style-type: none"> <li>● Unsatisfactory sample (saliva)</li> <li>● Sample delivery - <math>\geq 48</math> hours after sample collection</li> <li>● Storage temperature - <math>\geq 2-8^{\circ}\text{C}</math></li> </ul>	<ul style="list-style-type: none"> <li>● Low bacterial Load</li> <li>● normal flora overgrowth</li> <li>● normal flora overgrowth</li> </ul>
		Gastric lavage	<ul style="list-style-type: none"> <li>● No sodium bicarbonate added</li> <li>● Storage temperature - <math>\geq 2-8^{\circ}\text{C}</math></li> <li>● Sample delivery - <math>\geq 4</math> hours after sample collection</li> </ul>	<ul style="list-style-type: none"> <li>● acidic sample can inhibit bacterial growth</li> <li>● acidic sample can inhibit bacterial growth</li> <li>● acidic sample can inhibit bacterial growth</li> </ul>
	Line Probe Assay (LPA)	Sputum	<ul style="list-style-type: none"> <li>● AFB direct smear <math>&lt;1+</math></li> <li>● Sample delivery - <math>\geq 48</math> hours after sample collection</li> <li>● Storage temperature - <math>\geq 2-8^{\circ}\text{C}</math></li> </ul>	<ul style="list-style-type: none"> <li>● Low bacterial Load</li> <li>● normal flora overgrowth</li> <li>● normal flora overgrowth</li> </ul>
	PCR TB/NTM	Sputum	<ul style="list-style-type: none"> <li>● Sample delivery - <math>\geq 48</math> hours after sample collection</li> <li>● Storage temperature - <math>\geq 2-8^{\circ}\text{C}</math></li> <li>● Bloody sample</li> </ul>	<ul style="list-style-type: none"> <li>● normal flora overgrowth</li> <li>● normal flora overgrowth</li> <li>● Result invalid -RBC as an inhibitor</li> </ul>
	Interferon-gamma release assays (IGRA)	Whole blood	<ul style="list-style-type: none"> <li>● Insufficient sample – volume <math>&lt;0.8</math> ml or <math>&gt;1.2</math> ml</li> <li>● Sample hemolyzed</li> <li>● Sample lipemic</li> <li>● Improper sample shaking</li> <li>● Improper sample incubation</li> <li>● Improper sample centrifugation</li> </ul>	<ul style="list-style-type: none"> <li>● Result invalid</li> </ul>
	GeneXpert	Sputum	<ul style="list-style-type: none"> <li>● Unsatisfactory sample (saliva)</li> <li>● Sample delivery - <math>\geq 48</math> hours after sample collection</li> <li>● Storage temperature - <math>\geq 2-8^{\circ}\text{C}</math></li> <li>● Sample from patients who have received no antituberculosis therapy, or less than 3 days of therapy in the last 6 months.</li> <li>● Bloody sample – RBC as an inhibitor</li> </ul>	<ul style="list-style-type: none"> <li>● Low bacterial Load</li> <li>● normal flora overgrowth</li> <li>● normal flora overgrowth</li> <li>● Result interference (false result)</li> <li>● Result invalid -RBC as an inhibitor</li> </ul>
		Gastric lavage	<ul style="list-style-type: none"> <li>● Storage temperature - <math>\geq 2-8^{\circ}\text{C}</math></li> <li>● Sample delivery - <math>\geq 4</math> hours after sample collection</li> <li>● Sample from patients who have received no antituberculosis therapy, or less than 3 days of therapy in the last 6 months.</li> </ul>	<ul style="list-style-type: none"> <li>● acidic sample can inhibit bacterial grow</li> <li>● acidic sample can inhibit bacterial grow</li> <li>● Result interference (false result)</li> </ul>

**Note: Refer to Disease Section Test Handbook for sample collection and criteria**

UNIT	TEST	SAMPLE	FACTOR	CONSEQUENCES
Biokimia	TSH/FT4 Test Cholinesterase Test Liver Function Test Renal Profile Lipid Profile Miscellaneous Tests: -Calcium -Phosphorus -Fasting Blood Glucose -Random Blood Glucose	Serum /Plasma	<ul style="list-style-type: none"> <li>Inappropriate sample collection techniques for example drawing blood too quickly or leaving the tourniquet on for too long</li> <li>Improper order of draw</li> <li>Use of incorrect collection tubes</li> <li>Insufficient sample volume</li> <li>The collection site is not properly disinfected or if the collection tubes are contaminated.</li> <li>Inappropriate use of anticoagulants/ Improper mix of anticoagulant with sample or improper handling of blood</li> <li>Delay in transportation of samples</li> <li>Improper storage temperature of samples</li> <li>Improper &amp; delay of sample centrifugation/ separation</li> <li>Not fasting prior to sample collection of Fasting Blood Glucose Test &amp; Lipid Profile</li> </ul>	<ul style="list-style-type: none"> <li>Hemolysis (breakdown of red blood cells)</li> <li>Increase possibility of the carryover of additive contamination</li> <li>Inaccurate result/Hemolysis</li> <li>Inaccurate result</li> <li>Contamination of the sample</li> <li>Clotting</li> <li>Degradation of sample /Inaccurate result</li> <li>Degradation of sample /Inaccurate result</li> <li>Degradation of sample /Inaccurate result</li> <li>Inaccurate result</li> </ul>
	Adenosine Deaminase Test	Pleural fluid	<ul style="list-style-type: none"> <li>Incorrect sample collection technique (Aseptic technique during thoracentesis is required).</li> <li>Puncture Accidents- Accidental punctures during thoracentesis</li> <li>Insufficient sample volume</li> <li>Improper use of collection tube</li> <li>Delay in transportation of samples</li> <li>Improper storage temperature of samples</li> <li>Improper &amp; delay of sample centrifugation/ separation</li> </ul>	<ul style="list-style-type: none"> <li>Contamination of sample</li> <li>Can introduce blood into the sample</li> <li>Inaccurate result</li> <li>Inaccurate result</li> <li>Degradation of sample /Inaccurate result</li> <li>Degradation of sample /Inaccurate result</li> <li>Degradation of sample /Inaccurate result</li> </ul>
Bakteriologi	Culture test	Rectal swab/ Stool swab/ Throat swab/ Nasal swab/ Nasopharyngeal swab/ Pernal swab	<ul style="list-style-type: none"> <li>Improper storage temperature of samples (no ice)</li> <li>Delay in shipment (Sample delivery <math>\geq</math> 6 hours after sample collection)</li> <li>Wrong transport media</li> <li>Wrong swab</li> </ul>	<ul style="list-style-type: none"> <li>Overgrowth of Normal flora</li> <li>Low yield of pathogen &amp; Overgrowth of Normal flora</li> <li>No yield of pathogen</li> <li>No yield of pathogen</li> </ul>
	Leptospira MAT	Serum	<ul style="list-style-type: none"> <li>Improper &amp; delay of sample centrifugation/ separation</li> <li>Insufficient sample volume</li> </ul>	<ul style="list-style-type: none"> <li>Inaccurate result/ unable to see reaction with Leptospira serovar</li> <li>Unable to process sample</li> </ul>
	Leptospira PCR	Blood in EDTA	<ul style="list-style-type: none"> <li>Sent serum sample</li> <li>Insufficient sample volume</li> </ul>	<ul style="list-style-type: none"> <li>Leptospira DNA Not Detected</li> </ul>
	Syphilis Rechecking	Serum	<ul style="list-style-type: none"> <li>Improper &amp; delay of sample centrifugation/ separation</li> <li>Lysed serum</li> <li>Insufficient sample volume</li> </ul>	<ul style="list-style-type: none"> <li>Inaccurate result / unable to see agglutination</li> </ul>

**Note: Refer to Disease Section Test Handbook for sample collection and criteria**

UNIT	TEST	SAMPLE	FACTOR	CONSEQUENCES
Virology	Chikungunya/ Dengue/ Flavivirus/Zika Virus PCR	Serum /Blood	<ul style="list-style-type: none"> <li>Inadequate Sampel volume</li> <li>Sample not collected within specific days from onset of illness</li> <li>Storage temperature:&gt; 8°C</li> </ul>	<ul style="list-style-type: none"> <li><b>Inadequate sample size will result in template DNA</b> present in the reaction leading to <b>low yield of PCR product</b></li> <li>Decline in viral load</li> <li>Nucleic acids degradation</li> </ul>
	Enterovirus PCR –HFMD/ Conjunctivitis  Respiratory Virus PCR  (Influenza A H1N1, H3N2, Influenza B, Adenovirus, RSV, SARS- CoV-2, MERS- CoV)  Monkeypox Virus PCR	Mouth Ulcer Swab/ Vesicle Swab/ Rectal Swab/ Throat Swab/ Stool/ Pleural fluid / Eye swab (for Conjunctivitis) Nasopharyngeal Swab/ Nasal swab/ Throat Swab/ combined OPS/NPS  Lesion fluid swab /Scab or crust Tonsillar swab/ Nasopharyngeal swab	<ul style="list-style-type: none"> <li>Sample taken using other than flocced or dacron swab in VTM (volume 2-2.5ml) or sterile container.</li> <li>Sample not collected within specific days from onset of illness</li> <li>Storage temperature:&gt; 8°C</li> </ul>	<ul style="list-style-type: none"> <li>Flocced and Dacron swabs are non-inhibitory to PCR</li> <li>Decline in viral load</li> <li>Nucleic acids degradation</li> </ul>
	Rotavirus	Fresh stool	<ul style="list-style-type: none"> <li>Sample not collected within specific days from onset of illness</li> <li>Sent in unsuitable storage temperature</li> </ul>	<ul style="list-style-type: none"> <li>Decline in viral load</li> </ul>
	Hepatitis B surface antigen or antibody / Rubella /Measles IgM	Serum/Blood	<ul style="list-style-type: none"> <li>Sampel volume: &lt;1ml</li> <li>Lysed serum/Plasma</li> <li>Unsuitable storage temperature:8°C (&lt;24hours) OR &gt;-70°C (&gt;48 hours)</li> </ul>	<ul style="list-style-type: none"> <li>Inadequate sample for testing</li> <li>False positive results</li> <li>Degradation of sample</li> <li>Inaccurate result</li> </ul>
Parasite	PCR Malaria	Blood in EDTA	<ul style="list-style-type: none"> <li>Sample insufficient</li> <li>Send plasma or serum instead of EDTA</li> </ul>	<ul style="list-style-type: none"> <li>Inadequate sample for testing</li> <li>Inappropriate sample will lead to poor yield of DNA and lead to false negative result.</li> </ul>
	Slide BFMP	Slides	<ul style="list-style-type: none"> <li>Inappropriate packaging of slides</li> </ul>	<ul style="list-style-type: none"> <li>Potential of broken slides and unable to proceed with rechecking.</li> </ul>

**Note: Refer to Disease Section Test Handbook for sample collection and criteria**